



## High Risk Nitrile Examination Gloves



### Glove of Choice for First Responders, EMS/Police, Health Care/First Aid, Chemotherapy/Home Infusion, Compounding

- Meets **ASTM D6978** for Resistance to Permeation by **Chemotherapy Drugs**  
Permeation by **Fentanyl**
- Use for **USP800 Compliance**
- **High Risk** 8mil, 295mm/11.6" Extended Cuff
- Meets ASTM F1671 for Viral Penetration  
**0% Viral Penetration**  
**0% Alcohol Permeation**
- **Low Dermatitis** Patented Technology
- Quality Alternative to Latex Gloves

#### Tested for Use with **15** Chemotherapy Drugs



#### Superior Permeation Protection

- Meets **ASTM D6978** for Resistance to
  - Permeation by **Chemotherapy Drugs**
  - Permeation by **Fentanyl**

#### Minimum Breakthrough Detection Time in Minutes

##### Chemotherapy Drug & Concentration

> 240 minutes	Cisplatin 1.0 mg/ml, Cyclophosphamide (Cytoxan) 20 mg/ml, Cytarabine 100 mg/ml, Dacarbazine (DTIC) 10.0 mg/ml, Doxorubicin Hydrochloride 2.0 mg/ml, Etoposide 20.0 mg/ml, Fluorouracil 50.0 mg/ml, Ifosfamide 50.0 mg/ml, Methotrexate 25 mg/ml, Mitomycin C 0.5 mg/ml, Mitoxantrone 2.0 mg/ml, Paclitaxel (Taxol) 6.0 mg/ml, Vincristine Sulfate 1.0 mg/ml
---------------	---

##### Not Recommended for Use

Carmustine (BCNU) (3.3 mg/ml), Thiotepa (10.0 mg/ml)

##### Fentanyl Citrate & Concentration

> 240 minutes Fentanyl Citrate Injection (100mcg/2ml)

#### High Risk Nitrile Certification



Meets ASTM D6978



Use for USP800



Patented Technology\*  
**Low Derma**



#### UltiMates High Risk Nitrile Medical Exam Gloves - Product Numbers

##### 50ct - 8mil, 295mm/11.6" Extended Cuff

Small	<b>05631080</b>
Medium	<b>05632080</b>
Large	<b>05633080</b>
X-Large	<b>05634080</b>
XX-Large	<b>05635080</b>

NOTE: UltiMates glove packaging in English and Spanish.