



Find the Right Fit

CareMates Gloves should fit closely but should not be too tight. Prolonged use of gloves that are too tight-fitting may eventually contribute to health problems including joint pain, thumb fatigue, and carpal tunnel syndrome.



How to Find the Right Size

To determine the right size, align hand where thumb joins palm. Measure the widest part of your hand and line up with the sizes to the right.



NITRILE		VYTRILE		VINYL		LATEX	
100ct (100 gloves/box)		100ct (100 gloves/box)		100ct (100 gloves/box)		100ct (100 gloves/box)	
Small	3961265	Small	4763090	Small	3961083	Small	2533073
Medium	3961281	Medium	4090932	Medium	3961091	Medium	2533081
Large	3961299	Large	4565438	Large	3961133	Large	2533099
X-Large	4763157	X-Large	4740148	X-Large	(use Vytrile)	X-Large	(use Vytrile)
50ct (50 gloves/box)		50ct (50 gloves/box)		50ct (50 gloves/box)		50ct (50 gloves/box)	
Small	4273280	Small	3974425	Small	(use Vytrile)	Small	3757143
Medium	4273306	Medium	3974433	Medium	(use Vytrile)	Medium	3757150
Large	4273314	Large	3974466	Large	(use Vytrile)	Large	3757168



CardinalHealth

NOTE: To order, use the Cardinal CIN numbers listed above.